



NAME _____ PHONE NO _____
ACCOUNT # _____ P.O. _____ DATE _____
ADDRESS _____

**Note: This form is only for the ordering wheelchairs for the Manitoba Wheelchair Program.
A complete order form with all configurable options and accessories is available at www.invacare.ca**

BASE

TDXSP comes standard with: 300 lbs. user weight capacity, 5.8 MPH top speed, 4 pole motor package, MK6i SPJPLUS joystick, MK690 Controller, Quiet Stability Lock. Additional joystick options available.

- TDXSP TDX SP with MK6i™ Electronics
- TDXSPSEAT TDX Seat Only (For SMD Use Only) See Seat Selection
- TDXSPBASE FOR BASE ONLY, SEE TDXSP-CG Form #06-490C

TRANSPORT TIE-DOWN*

TRBKTS includes four factory installed wheelchair transport brackets. TRBKTS has not been crash-tested in accordance with WC19. Use these transport brackets only to secure an unoccupied wheelchair during transport. As of this date, the Department of Transportation has not approved any tie-down systems for transportation of a user while in a wheelchair, in a moving vehicle of any type. It is Invacare's position that users of wheelchairs should be transferred into appropriate seating in vehicles for transportation and use be made of the restraints made available by the auto industry. Invacare cannot and does not recommend any wheelchair transportation systems.

* Transport Tie Down option is not retrofittable to existing models and is not field serviceable.

TRANSPORT TIE DOWN

- TRBKTS Wheelchair Transport Brackets

USER WEIGHT LIMITS

- U300 User Weight <300 lbs.
Package includes 4-pole motors and 5.8 MPH top speed with <300 lbs. user weight limit.

Note: For G-Trac Module, go to Controller area under MK6i Selection

SEAT TO FLOOR HEIGHTS

Due to variation in options, the listed seat-to-floor heights have a tolerance of +/- .5"

- LOW Seat-to-Floor Low (16.5" at 0°)
- MED Seat-to-Floor Medium (18.5" at 0°)
- TALL Seat-to-Floor Tall (20.5" at 0°)

WHEEL OPTION (PAIR) *Note: SMD will not cover the cost of servicing pneumatic tires

- 1430-3T 14" X 3" Pneumatic Tires
- 1431-3T 14" X 3" Tire with Foam Filled Inserts

FRAME FINISH

- 24P Wet Black

SEAT

SEAT SELECTION

- ADJASBA Adjustable Rehab Seat with Adjustable Seat Tilt & Back Angle

SEAT WIDTH RANGE

- 1620W 16" to 20" Seat Width Range
- 2024W 20" to 24" Seat Width Range

FOOTNOTES:
1 - Option to purchase.

SEAT WIDTH

Contoura available in 16", 18", 20" and 22" seat widths.

- W16 Seat Width 16"
- W17 Seat Width 17"
- W18 Seat Width 18"
- W19 Seat Width 19"
- W20 Seat Width 20"
- W21 Seat Width 21"
- W22 Seat Width 22"

SEAT DEPTH RANGE

- 1619D 16" to 19" Seat Depth Range
- 1922D 19" to 22" Seat Depth Range

SEAT DEPTH SETTING

- D16 Seat Depth 16"
- D17 Seat Depth 17"
- D18 Seat Depth 18"
- D19 Seat Depth 19"
- D20 Seat Depth 20"
- D21 Seat Depth 21"
- D22 Seat Depth 22"

ARM TYPE RIGHT

- 29FR Desk Length Adjustable Height Right
- 19FR Full Length Adjustable Height Right

ARM TYPE LEFT

- 19FL Full Length Adjustable Height Left
- 29FL Desk Length Adjustable Height Left

ASBA SEAT BACK TYPE

- SECT Sectional Rear Uprights
- L8661 Omit Rear Uprights

ASBA BACK ANGLE DEGREE

- 90BA 90°
- 95BA 95°

ASBA SEAT BACK HEIGHT

- BH16 Back Height 16"
- BH17 Back Height 17"
- BH18 Back Height 18"
- BH19 Back Height 19"
- BH20 Back Height 20"
- BH21 Back Height 21"
- BH22 Back Height 22"
- BH23 Back Height 23"
- BH24 Back Height 24"

CONTOURA BACK HEIGHT

- TSBH22 Back Height 25"
- TSBH25 Back Height 25"
- TSBH28 Back Height 28"

SEAT POSITIONING STRAP

- 1515 Auto Style Seat Positioning Strap

CHEST POSITIONING STRAP

- 7321BK Hook & Loop Chest Positioning Strap

FOOTNOTES:

- 1 - Option to purchase.
- 2 - Justification required.

ELECTRONICS

MK6i™ CONTROLLERS

- MK690 90 AMP Controller
- GTRAC G-Trac™ Electronic Gyroscope Module \$988 (1)

G-Trac provides Improved Capabilities and enhanced Traction. Does not require Special Motors. Not available with SPJPLUS joysticks.

STANDARD DRIVER CONTROL

- SPJPLUS SPJPLUS Single Drive Joystick

MULTIPLE DRIVE CONTROL OPTIONS

- PWH Harness required for Expandable Electronics
Harness is required to interface joystick or MK6i display with expandable controllers. Option is located under Electronics Accessories in configurator.
- CMPJM6U MPJ+ Multiple Drive Joystick with Colour Display (2)
- ASLRDYM6 MK6 Display Only (No Driver Control) (2)
Display Only, No Driver control. Allows choice from order form.

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS

- 1558MM6 Compact Joystick with Display (2)

ALTERNATIVE DIGITAL CONTROL SYSTEMS

Must also use either a Multiple Drive Joystick or MK6i Display.

- SNPM6 Sip-N-Puff Digital Interface Box (2)

SIP N PUFF KIT

- PKG32666 Therafin Sip-N-Puff Breath Tube Kit (2)

ADDITIONAL CONTROL CHOICES

- 1558M6 Compact Joystick - Less Standard Mount \$937 (1)
- PACM6 Proportional Attendant Control with Quick Release Mount \$1,024 (1)

JOYSTICK MOUNTING HARDWARE

- QLAM6 Quad Link Retractable Joystick Mount CMPJ+ \$358 (1)
QLAM6 and QLA come standard with Stealth Adj. Height Joystick Mount
- QLA Quad Link Retractable Joystick Mount SPJPLUS \$358 (1)
QLAM6 and QLA come standard with Stealth Adj. Height Joystick Mount
- LEFT Left Hand Mounted Joystick

ELECTRONIC ACCESSORIES

- 1813M6 MK6i Programmer with Professional Memory Card
- 1813SD Professional Memory Card
- A24VPS 24V MK6i Auxiliary Power Source \$252 (1)
- AUX12M6 Auxiliary Module For 1 & 2 \$706 (1,2)
- AUX34M6 Auxiliary Module For 3 & 4 \$706 (1,2)
- MPY Mono Port "Y" Cable (2)
For use with MPJ+ or MK6i Display to provide additional switch sites.
- 2PB Dual Push Button
Option to replace Mono Port "Y" cable when requiring two switch functions of the mode switch port on expandable controller systems. Compatible with CMPJM6U. Not compatible with SPJPLUS.
- PBMT Push Button Mounting Kit

Memory Card - BASIC

- Standard with ALL Expandable Control Systems (CMPJ+, MK6i Display).
- Used to back up / restore programmed settings for only one chair.
- May be given to user or kept with provider's files for safe keeping.
- Does not contain advanced diagnostics / help screens / file structure.
- Not compatible with SPJ+ Joysticks.

Memory Card - PROFESSIONAL

- Standard with ALL MK6i Programmers.
- Available also with USB Card Reader on Order Form.
- Contains advanced diagnostics, help screens, software upgrades, file storage / retrieval
- Used to back up multiple chairs and programming settings.
- Intended for use by qualified providers only.
- Not compatible with SPJ+ Joysticks.

MISCELLANEOUS

BATTERY TRAY TYPE

22TRY Base Width 24"; 24TRY Base Width 25.5"

- 22TRY 22 NF Style Tray
- 24TRY 24 NF Style Tray

INSTALLATION OF BATTERIES

- 22NFBATTERY 22 NF Battery (On Chair)
- 24BATTERY 24 Gel Battery (On Chair) (2)

BATTERY CHARGER OPTIONS

- 110CHARGER 110 Volt Battery Charger

BACK ACCESSORIES

- L8601 Omit Back Upholstery
- RBP Reinforced Back Upholstery

FOOTNOTES:
1 - Option to purchase.

MISCELLANEOUS ACCESSORIES

- 1303 Wheel Locks \$93 (1)

FRONT RIGGING

All Front rigging options include a standard footplate with the exception of the 70NHD.

HANGER BRACKET TYPE

- ASB Center Pivot Style Bracket

CENTRE PIVOT STYLE

- 70NHD 70° Non-Tapered S/A Footrest 13.25"-17" (Steel)
Comes standard with 1651 footplates. With seat widths of 21" to 22" and depth of 16", must pick footplate 1350. The AT5543 is not available.
- 70STAPER 70° Tapered Footrests 13"-17" (Aluminum)
Standard with impact guards. Not available with widths 19" - 24" and 16" deep seat. Comes with composite footplates.
- AT5044 Elevating Painted Legrests - Wedge Style (2)
Not available with 16" wide by 16" to 18" deep or 17" wide by 16" deep and low seat-to-floor.

CENTRE MOUNT

- CTMT-F Centre Mount Adjustable Knee Angle 12 1/4"-15"
Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 12.5" - 15"
- CTMT-FJ Junior Centre Mount Adjustable Knee Angle 9 1/4"-12"
Allows 60 to 90 degree static knee angle positioning with a drop back lever for easier transferring. Extension range 9.25" - 12"
- SM Small Footplates (5.25" x 7.5")
- MED Medium Footplates (5.50" x 9.5")
- LRG Large Footplates (5.90" x 11.5")
- CALF Calf Pads For Center Mount \$268 (1)

FOOTPLATE OPTIONS (PAIR)

- 1651 Composite Footplates
Not available with AT5044 and 17" wide by 16" deep and low seat-to-floor.
- 1350 Extra Large Aluminum Footplates (2)
Not available with AT5044 and 16" wide by 16"-18" or 17" wide by 16" deep and low seat-to-floor.
- AT5543 Adjustable Angle Flip Up Footplates (2)

FOOTREST/LEGREST ACCESSORIES

- 1337 Calf Strap
- HLOOP Heel Loops
- IMPACTG Impact Guards
- LPT2 Longer Pivot Slide Tube For 1350
- ALPT Longer Pivot And Slide Tube
- ALPT4 Longer Pivot and Slide Tube 4"

FOOTNOTES:
1 - Option to purchase.

TDX SP Specifications		
	Adult (Adjust)	
Seat Width Range (adjustable in 1" increments)	16" - 22"	
Seat Depth Range	16" - 22"	
Back Height Range	16" - 24"	
Back Angle Range	85° - 105° in 5° increments	
Seat-to-Floor (approx.)		
Standard	16 1/2" at 0°	
Medium	18 1/2" at 0°	
Tall	20 1/2" at 0° at 5° add 1/2"	
Overall Height (with 20" back height and 0° seat angle)		
Standard	35 1/2"	
Medium	37 1/2"	
Tall	39 1/2"	
Overall Length	35 1/4" without front rigging	43" with 90° front rigging
Weight		
4-Pole Motor	without 22NF Batteries	166 lbs.
	with two 22NF Batteries	264 lbs.
Drive Axle	Non-Adjustable	
Drive Wheels/Tires	(Foam-Filled or Pneumatic) Standard	14" x 3"

Weight Limit 300 Pounds
Invacare recommends all wheelchairs be fitted and programmed to the consumer under the guidance of a healthcare professional, knowledgeable of the particular limitations of the product.

Note: All Weight limits are inclusive of add-on equipment and seating systems. All specifications and dimensions are approximate.

Batteries
22NF batteries with terminal configuration (positive on the left and negative on the right) MUST be used. 22NF batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection.

GP24 batteries with terminal configuration (negative on the left and positive on the right) MUST be used. GP24 batteries that have the reverse terminal configuration MUST not be used. Terminals MUST have a cross hole for proper battery connection. See Owner's Manual, part number 1114809. These recommendations MUST be followed otherwise injury and damage may occur.

TDX SP Extended Warranty Information
Electrical components, motors - 2 years
Base frame - 5 years
Batteries - 6 months
*Warranty is valid to the original purchaser who purchases this product within Canada.
*For more product warranty information, please refer to the product owners' manual.

Invacare Canada reserves the right to change prices and specifications without notice.
HEAD OFFICE - 570 Matheson Blvd E., Unit #8, Mississauga, ON L4Z 4G4 Phone: (905) 890-8300 (800) 668-5324 (Canada) Fax: (905) 890-5244 (800) 668-5478 (Canada)
Regional Centre - 16769 Boul. Hymus, Kirkland, P.Q. H9H 3L4 Phones: (514) 630-6080 (800) 363-0137(P.Q.) Fax: (514) 630-6120
Invacare®, and Yes, you can., are trademarks of Invacare Corporation.
WWW.INVACARE.CA

NOTE: All specifications and dimensions are approximate
For the most current pricing information visit www.invacare.ca